

Technical Datasheet

Article - No:	TDS-RT050TFX-L	Revision: A0
Remarks	single use only	
Specification		
Product name	RT050TFX-L, 50µL Tecan robotic conductive filter tips, low retention	
Production conditions	Class 100,000 clean room production	
	No release agent is used	
	Motor molding	
Classification	CE/IVD	
Dimensions		
Colour	Black	
Nominal Capacity	1-50µL	
Material Characteristics	polypropylene (PP), premium surface filter: high density polyethylene (PE)	
Product Characteristics	biologically inert	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogen	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
Match the model	free of ATP	
	Genesis Freedom [®]	
	Freedom EVO [®]	
Shelf - Life	Miniprep with LiHa [®]	
	36 months at controlled storage conditions	
Packaging	pcs per rack: 96	
	rack per case: 30	
Storage	Storage Temperature: 2°C-40°C, humidity 40%-80%	
	Keep clean, ventilated and out of direct sunlight	
Approved by PakGent quality department		