

Technical Datasheet

Article - No:	TDS-UFPT-1000R	Revision: A0
Remarks	single use only	
Specification		
Product name	UFPT-1000R, 1000µL Pipette tips	
Production conditions	Class 100,000 clean room production	
	No release agent is used	
	Motor molding	
Classification	CE/IVD	
Dimensions	<p>Technical drawing of a pipette tip. The drawing shows a tapered cylindrical shape. The total length is indicated as 78⁺⁰. The outer diameter at the wider top end is 9.7. The outer diameter at the narrower bottom end is 1.1±0.05. The inner diameter at the bottom end is 0.72±0.05. There are also two diameters indicated at the top end: 8.5 and 9.7.</p>	
Colour	Blue	
Graduation	moulded rings	
Nominal Capacity	100-1000µL	
Material Characteristics	polypropylene (PP), premium surface	
Product Characteristics	biologically inert	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogen	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
Shelf - Life	36 months at controlled storage conditions	
Packaging	pcs per rack: 96	
	rack per case: 50	
Storage	Storage Temperature: 2°C-40°C, humidity 40%-80%	
	Keep clean, ventilated and out of direct sunlight	
Approved by PakGent quality department		