

Technical Datasheet

Article - No:	TDS-RT1200R	Revision: A0
Remarks	single use only	
Specification		
Product name	RT1200R, 1200µL Rainin LTS Pipette tips	
Production conditions	Class 100,000 clean room production	
	No release agent is used	
	Motor molding	
Classification	CE/IVD	
Dimensions	<p>The diagram shows a pipette tip with a total length of 86.6. The top left diameter is 8.1, and the top right diameter is 8.8. The volume is 20µL. The diameter at the tip is 1.34±0.05, and the diameter at the base is 0.95±0.05.</p>	
Colour	Clear	
Graduation	moulded rings	
Nominal Capacity	100-1200µL	
Material Characteristics	polypropylene (PP), premium surface	
Product Characteristics	biologically inert	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogen	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
Shelf - Life	36 months at controlled storage conditions	
Packaging	pcs per box: 96	
	box per case: 50	
Storage	Storage Temperature: 2°C-40°C, humidity 40%-80%	
	Keep clean, ventilated and out of direct sunlight	
Approved by PakGent quality department		