

Technical Datasheet

Article - No:	TDS-RT1000T	Revision: A0
Remarks	single use only	
Specification		
Product name	RT1000T, 1000µL Tecan robotic pipette tips	
Production conditions	Class 100,000 clean room production	
	No release agent is used	
	Motor molding	
Classification	CE/IVD	
Dimensions	<p>The drawing shows a tapered pipette tip with a cylindrical base. Dimension 1 is the total length (95.60 ± 0.25). Dimension 2 is the diameter of the cylindrical base (5.00 ± 0.0). Dimension 3 is the diameter of the narrow tip section (7.90 ± 0.0). Dimension 4 is the thickness of the base (1.00 ± 0.05).</p>	
Colour	Clear	
Nominal Capacity	100-1000µL	
Material Characteristics	polypropylene (PP), premium surface	
Product Characteristics	biologically inert	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogen	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
	free of ATP	
Match the model	Genesis Freedom [®]	
	Freedom EVO [®]	
	Miniprep with LiHa [®]	
Shelf - Life	36 months at controlled storage conditions	
Packaging	pcs per rack: 96	
	rack per case: 60	
Storage	Storage Temperature: 2°C-40°C, humidity 40%-80%	
	Keep clean, ventilated and out of direct sunlight	
Approved by PakGent quality department		