

Technical Datasheet

| | | |
|---|---|--------------|
| Article - No: | TDS-RT300-REF | Revision: A0 |
| Remarks | single use only | |
| Specification | | |
| Product name | RT300-REF,300µL Pipette tips,Refull pack | |
| Production conditions | Class 100,000 clean room production | |
| | No release agent is used | |
| | Motor molding | |
| Classification | CE/IVD | |
| Dimensions | | |
| Colour | Clear | |
| Graduation | moulded rings | |
| Nominal Capacity | 1-300µL | |
| Sterile | Yes | |
| Material Characteristics | polypropylene(PP),premium surface pack:Polyethylene terephthalate(PET) | |
| Product Characteristics | biologically inert | |
| | no heavy metals according to EC/94/62 | |
| Purity Grade | free of RNase and DNase | |
| | free of pyrogen | |
| | free of PCR inhibitors | |
| | free of human DNA | |
| | free of endotoxins | |
| | free of ATP | |
| Shelf - Life | 36 months at controlled storage conditions | |
| Packaging | pcs per rack: 96 | |
| | rack per pack: 8 | |
| | pack per case: 10 | |
| Storage | Storage Temperature:2°C-40°C,humidity40%-80% | |
| | Keep clean,ventilated and out of direct sunlight | |
| Approved by PakGent quality department | | |