


Technical Datasheet

Article - No:	TDS-STSS-510CB	Revision: A0
Remarks	single use only	
Specification		
Product name	STSS-510CB, 1.5ml Screw Cap Tubes	
Production conditions	Class 100,000 clean room production	
	No release agent is used	
Classification	CE/IVDD	
Dimensions		
Colour	Clear tube, Blue cap	
Nominal Capacity	1.5ml	
Material Characteristics	polypropylene (PP)	
Performance parameters	Storage Temperature: -80°C~121°C	
Product Characteristics	biologically inert	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogen	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
Shelf - Life	36 months at controlled storage conditions	
Packaging	pcs per bag: 500	
	bag per case: 10	
Storage	dry condition: 40%-60% rel. humidity	
Approved by PakGent quality department		