

Technical Datasheet

Article - No:	TDS-RT050TX	Revision: A0
Remarks	single use only	
Specification		
Product name	RT050TX,50µL Tecan robotic conductive tips	
Production conditions	Class 100,000 clean room production	
	No release agent is used	
	Motor molding	
Classification	CE/IVD	
Dimensions		
Colour	Black	
Nominal Capacity	1-50µL	
Material Characteristics	polypropylene (PP),premium surface	
Product Characteristics	biologically inert	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogen	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
	free of ATP	
Match the model	Genesis Freedom [®]	
	Freedom EVO [®]	
	Miniprep with LiHa [®]	
Shelf - Life	36 months at controlled storage conditions	
Packaging	pcs per rack: 96	
	rack per case: 30	
Storage	Storage Temperature:2°C-40°C,humidity40%-80%	
	Keep clean,ventilated and out of direct sunlight	
Approved by PakGent quality department		