

Technical Datasheet

Article - No:	TDS-UFPT-10EB-L	Revision: A0
Remarks	single use only	
Specification		
Product name	UFPT-10EB-L, 10 μ L Pipette tips, low retention	
Production conditions	Class 100,000 clean room production	
	No release agent is used	
	Motor molding	
Classification	CE/IVD	
Dimensions	<p style="text-align: center;"> $\phi 4.3$ $\phi 5.9$ 46.0 $10\mu\text{L}$ $\phi 0.84 \pm 0.05$ $\phi 0.46 \pm 0.05$ </p>	
Colour	Clear	
Graduation	moulded rings	
Nominal Capacity	0.1-20 μ L	
Material Characteristics	polypropylene (PP), premium surface	
Product Characteristics	biologically inert	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogen	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
Shelf - Life	36 months at controlled storage conditions	
Packaging	pcs per bag: 1000	
	bag per case: 40	
Storage	Storage Temperature: 2°C-40°C, humidity 40%-80%	
	Keep clean, ventilated and out of direct sunlight	
Approved by PakGent quality department		