

Technical Datasheet

Article - No:	TDS-CL-D100	Revision: A0
Remarks	single use only	
Specification		
Product name	CL-D100,100mm,Cell Culture Dish	
Production conditions	Class 10,000 clean room production	
	No release agent is used	
Classification	CE/IVD	
Dimensions	<p>The technical drawings include: <ul style="list-style-type: none"> TOP View: A circular dish with a diameter of 92.2 mm. Side View: Shows a total height of 21.49 mm. The inner well has a diameter of 88.58 mm and a depth of 87.35 mm. Bottom View: Shows the bottom of the dish with a diameter of 94.93 mm. Detail View: A separate view of the bottom flange. </p>	
Colour	Clear	
Nominal Capacity	100mm	
Material Characteristics	polystyrene (PS), premium surface	
Sterile	Yes	
Product Characteristics	biologically inert	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogen	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
Shelf - Life	36 months at controlled storage conditions	
	pcs per bag: 10 bag per case: 30	
Storage	Storage Temperature:20°C-25°C,humidity40%-60%	
	Keep clean,ventilated and out of direct sunlight	
Approved by PakGent quality department		