



This is to certify that units of the lot number below were tested and found to comply with the specifications of the grade listed. Products are tested to meet the analytical requirements of the noted grade. The following information is the actual analytical results obtained, for reference.

Part Number:	PD150		
Lot.No.:	202505090101		
	Size: <input type="checkbox"/> 35mm <input type="checkbox"/> 60mm <input type="checkbox"/> 90mm <input type="checkbox"/> 100mm <input checked="" type="checkbox"/> 150mm		
Specification:	Surface treatment: <input type="checkbox"/> Yes <input checked="" type="checkbox"/> No		
Product Name:	Petri Dish		
Package:	10pcs/bag, 16bags/case		
Shelf Life:	3 years		
Manufacturer:	PakGent Bioscience (Suzhou) Co., Ltd.		
Inspection Date:	15th,May,2025	Report Date	15th,May,2025
Inspection items and result:			
Cosmetic:	No black spots, no cuts,clear and transparent		
Package:	Accurate quantity and label, inner package without damage		
Dimension:	Planeness Lid $\leq 0.5\text{mm}$, Dish $\leq 0.5\text{mm}$		
Function test:	N/A		
Endotoxin:	$<0.03\text{EU/ml}$		
RNase/DNase:	DNase $\leq 4 \times 10^{-6}$ Kunitz units/ μl ; RNase $\leq 10^{-9}$ Kunitz units/ μl		
PCR Inhibitors:	Not detectable		
Human DNA:	$<30\text{pg}$		
Sterile	EO sterile		
Conclusion: Comply with Specification			

This signature indicates that the above test has met all quality specification and has been reviewed by a quality representative.

Title: QA Manager

Signature: Alvin Cao