


## Technical Datasheet

<b>Article - No:</b>	TDS-DWPS20U	Revision: A0																				
<b>Remarks</b>	single use only																					
Specification																						
<b>Product name</b>	DWPS20U,2.2ml deep well plate,Square well U bottom																					
<b>Production conditions</b>	Class 100,000 clean room production																					
	No release agent is used																					
	Motor molding																					
<b>Classification</b>	CE/IVD																					
<b>Dimensions</b>	 <table border="1" style="margin-left: auto; margin-right: auto; border-collapse: collapse;"> <thead> <tr> <th>Part No.</th> <th>Description</th> <th>A ±0.2</th> <th>B ±0.2</th> <th>C ±0.2</th> <th>D ±0.1</th> <th>E ±0.3</th> <th>F ±0.05</th> <th>G ±0.2</th> <th>H +0.4/-0.2</th> </tr> </thead> <tbody> <tr> <td>DWPS20U</td> <td>2.2ml square well plate, round well U bottom</td> <td>127.30</td> <td>85.30</td> <td>41.60</td> <td>0.60</td> <td>2.95</td> <td>8.20</td> <td>123.13</td> <td>80.93</td> </tr> </tbody> </table>		Part No.	Description	A ±0.2	B ±0.2	C ±0.2	D ±0.1	E ±0.3	F ±0.05	G ±0.2	H +0.4/-0.2	DWPS20U	2.2ml square well plate, round well U bottom	127.30	85.30	41.60	0.60	2.95	8.20	123.13	80.93
Part No.	Description	A ±0.2	B ±0.2	C ±0.2	D ±0.1	E ±0.3	F ±0.05	G ±0.2	H +0.4/-0.2													
DWPS20U	2.2ml square well plate, round well U bottom	127.30	85.30	41.60	0.60	2.95	8.20	123.13	80.93													
<b>Colour</b>	Clear																					
<b>Bottle Shape</b>	U type																					
<b>Nominal Capacity</b>	2.2ml																					
<b>Material Characteristics</b>	polypropylene (PP),premium surface																					
<b>Product Characteristics</b>	biologically inert																					
	no heavy metals according to EC/94/62																					
<b>Purity Grade</b>	free of RNase and DNase																					
	free of endotoxins																					
<b>Shelf - Life</b>	36 months at controlled storage conditions																					
<b>Packaging</b>	pcs per bag: 10																					
	bag per case: 10																					
<b>Storage</b>	Storage Temperature:2°C-40°C,humidity40%-80%																					
	Keep clean,ventilated and out of direct sunlight																					
<b>Approved by PakGent quality department</b>																						