

Technical Datasheet

Article - No:	TDS-CL-D035	Revision: A0
Remarks	single use only	
Specification		
Product name	CL-D035,35mm,Cell Culture Dish	
Production conditions	Class 10,000 clean room production	
	No release agent is used	
Classification	CE/IVD	
Dimensions	<p>The drawing shows three views of the cell culture dish: <ul style="list-style-type: none"> TOP: A circular view of the dish. Side View: A cross-section showing the dish's profile with dimensions: an outer diameter of $\phi 37.95$, an inner diameter of $\phi 36.82$, and a height of 12.38. Bottom: A view of the dish's base with a diameter of $\phi 38.79$. </p>	
Colour	Clear	
Nominal Capacity	35mm	
Material Characteristics	polystyrene (PS), premium surface	
Sterile	Yes	
Product Characteristics	biologically inert	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogen	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
Shelf - Life	36 months at controlled storage conditions	
	pcs per bag: 20 bag per case: 25	
Storage	Storage Temperature:20°C-25°C,humidity40%-60%	
	Keep clean,ventilated and out of direct sunlight	
Approved by PakGent quality department		