

Technical Datasheet

| | | |
|---|--|--------------|
| Article - No: | TDS-EFBC-125P | Revision: A0 |
| Remarks | single use only | |
| Specification | | |
| Product name | EFBC-125P, 125ml, Erlenmeyer Flask, Plug Cap, Baffled bottom | |
| Production conditions | Class 100,000 clean room production | |
| | No release agent is used | |
| Classification | CE/IVD | |
| Dimensions | <p style="font-size: small;">螺纹为右旋单头1.5圈, 螺距4.23mm</p> | |
| Colour | Clear | |
| Graduation | moulded rings | |
| Nominal Capacity | 125ml | |
| Material Characteristics | bottle: Polycarbonate(PC) | |
| | cap: High Density Polyethylene(HDPE) | |
| Sterile | Yes | |
| Product Characteristics | biologically inert | |
| | no heavy metals according to EC/94/62 | |
| Purity Grade | free of RNase and DNase | |
| | free of pyrogen | |
| | free of PCR inhibitors | |
| | free of human DNA | |
| | free of endotoxins | |
| | free of ATP | |
| Shelf - Life | 36 months at controlled storage conditions | |
| Packaging | pcs per case: 24 | |
| Storage | Storage Temperature: 20°C-25°C, humidity 40%-60% | |
| | Keep clean, ventilated and out of direct sunlight | |
| Approved by PakGent quality department | | |