


Technical Datasheet

Article - No:	TDS-DWPS09V	Revision: A0
Remarks	single use only	
Specification		
Product name	DWPS09V,0.9 ml 96 well plate,square well V bottom	
Production conditions	Class 100,000 clean room production	
	No release agent is used	
	Motor molding	
Classification	CE/IVD	
Dimensions		
Colour	Clear	
Bottle Shape	V type	
Nominal Capacity	0.9ml	
Material Characteristics	polypropylene (PP), premium surface	
Product Characteristics	biologically inert	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of endotoxins	
Shelf - Life	36 months at controlled storage conditions	
Packaging	pcs per bag: 10	
	bag per case: 10	
Storage	Storage Temperature:2°C-40°C,humidity40%-80%	
	Keep clean,ventilated and out of direct sunlight	
Approved by PakGent quality department		