


Petri Dish


Petri dishes are one of the most important products in bacteriological laboratories. Manufactured from high-grade polystyrene, the petri dishes have exceptional optical clarity for microscopic analysis as well as heat resistance for use with hot agar. To enable parallel testing of a single sample using different agars, divided dishes with two or three or four compartments or quadrant grid are available.




| Features | Description |
|-----------------------------|---|
| Material | High-grade polystyrene |
| Compatibility | Compliant with automatic plate pourers |
| Specification | 35 mm, 60 mm, 90 mm, 100 mm, 150 mm |
| Appearance | Divided dishes with two or three or four compartments, quadrant grid, flat, are available |
| Sterilization method | EO |



Transparent and clean
easy to observe under a microscope



Stable stacking



Double-layer protective packaging for
easy use and transportation



| Part NO. | Description | Pack |
|----------|--|---------------------------|
| PD90 | Petri Dish, 90*15mm, PS, whithout external grip | 20 pcs/bag, 25 bags/case |
| PD60 | Petri Dish, 60*15mm, PS, whithout external grip | 20 pcs/bag, 50 bags/case |
| PD35 | Petri Dish, 35*15mm, PS, whithout external grip | 10 pcs/bag, 230 bags/case |
| PD90-2 | Petri Dish, 90*15mm, 2 compartments, PS, whithout external grip | 20 pcs/bag, 25 bags/case |
| PD90-3 | Petri Dish, 90*15mm, 3 compartments, PS, whithout external grip | 20 pcs/bag, 25 bags/case |
| PD90-4 | Petri Dish, 90*15mm, 4 compartments, PS, whithout external grip | 20 pcs/bag, 25 bags/case |
| PD60-RGC | Rodac Dish, quadrant grid, 60mm, PS, concave, whithout external grip | 20 pcs/bag, 50 bags/case |
| PD60-RGF | Rodac Dish, quadrant grid, 60mm, PS, flat, whithout external grip | 20 pcs/bag, 50 bags/case |
| PD100 | Petri Dish, 100*15mm, PS | 10 pcs/bag, 32 bags/case |
| PD150 | Petri Dish, 150*15mm, PS | 10 pcs/bag, 16 bags/case |

PakGent®

Save every drop of precious reagent for you.

PakGent Bioscience (Suzhou) Co., Ltd.
 Tel: +86 (0) 512 6670 0017
 Address: Floor 1-3, Building 33, Alloy Industrial Park, Wuzhong District, Suzhou, China.

PakGent Bioscience (Changshu) Co., Ltd.
 Tel: +86 (0) 512 5298 8330
 Address: No.9 Lecheng Road, Bixi Street, Changshu, Jiangsu Province, China.

PakGent Bioscience (Shanghai) Co., Ltd.
 Tel: +86 (0) 512 6670 0017
 Address: Building C, No 888, Lake west road, Lin gang section, Free Trade Zone, Shanghai, China.