

Technical Datasheet

Article - No:	TDS-UFPT-REF-10ET	Revision: A0
Remarks	single use only	
Specification		
Product name	UFPT-REF-10ET, 10µL Extra long Pipette tips, Refill pack	
Production conditions	Class 100,000 clean room production	
	No release agent is used	
	Motor molding	
Classification	CE/IVD	
Dimensions	<p style="text-align: center;"> 46 mm $\phi 3.9$ $\phi 6.4$ $10\mu\text{L}$ $2.5\mu\text{L}$ $\phi 0.76 \pm 0.05$ $\phi 0.36 \pm 0.05$ </p>	
Colour	Clear	
Graduation	moulded rings	
Nominal Capacity	0.1-20µL	
Sterile	Yes	
Material Characteristics	polypropylene(PP), premium surface	
	pack: Polyethylene terephthalate(PET)	
Product Characteristics	biologically inert	
	no heavy metals according to EC/94/62	
	free of RNase and DNase	
Purity Grade	free of pyrogen	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
	free of ATP	
	free of ATP	
Shelf - Life	36 months at controlled storage conditions	
Packaging	pcs per rack: 96	
	rack per pack: 9	
	pack per case: 10	
Storage	Storage Temperature: 2°C-40°C, humidity 40%-80%	
	Keep clean, ventilated and out of direct sunlight	
Approved by PakGent quality department		