


Technical Datasheet

Article - No:	TDS-PT-10000RF-E	Revision: A0
Remarks	single use only	
Specification		
Product name	PT-10000RF-E, 10mL filter tips	
Production conditions	Class 100,000 clean room production	
	No release agent is used	
	Motor molding	
Classification	CE/IVD	
Dimensions		
Colour	Clear	
Graduation	moulded rings	
Nominal Capacity	10mL	
Material Characteristics	polypropylene(PP), premium surface	
	filter: high density polyethylene(PE)	
Product Characteristics	biologically inert	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogen	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
	free of ATP	
Match the model	Eppendorf	
	Biohit	
	DLAB	
Shelf - Life	36 months at controlled storage conditions	
Packaging	pcs per rack: 24	
	rack per case: 40	
Storage	Storage Temperature: 2°C-40°C, humidity 40%-80%	
	Keep clean, ventilated and out of direct sunlight	
Approved by PakGent quality department		