

Technical Datasheet

Article - No:	TDS-DWPC96	Revision: A0																																																							
Remarks	single use only																																																								
Specification																																																									
Product name	DWPC96, KingFisher Plate Comb																																																								
Production conditions	Class 100,000 clean room production																																																								
	No release agent is used																																																								
	Motor molding																																																								
Classification	CE/IVD																																																								
Dimensions																																																									
	<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th colspan="2">Plate dimension 板尺寸[mm]</th> <th colspan="4">Well dimension 孔尺寸[mm]</th> </tr> <tr> <th colspan="2">Length 长度</th> <th colspan="2">Width 宽度</th> <th colspan="2">Height 高度</th> <th colspan="2">Depth 深度</th> </tr> <tr> <th>L</th> <th>L1</th> <th>3</th> <th>H</th> <th>H1</th> <th>h</th> <th>h1</th> <th></th> </tr> </thead> <tbody> <tr> <td>127.30</td> <td>110.00</td> <td>85.30</td> <td>45.00</td> <td>48.00</td> <td>43.00</td> <td>43.90</td> <td></td> </tr> <tr> <th colspan="2">Bottom 底部[mm]</th> <th colspan="4">Position of well center 孔位置[mm]</th> <th rowspan="3">Theoretical Weight[g]</th> </tr> <tr> <th>ID 内径</th> <th>OD 外径</th> <th></th> <th></th> <th></th> <th></th> </tr> <tr> <td>6</td> <td>D</td> <td>P1</td> <td>P2</td> <td></td> <td></td> </tr> <tr> <td>4.30</td> <td>5.60</td> <td>9.00</td> <td>9.00</td> <td></td> <td></td> <td></td> </tr> </tbody> </table>		Plate dimension 板尺寸[mm]		Well dimension 孔尺寸[mm]				Length 长度		Width 宽度		Height 高度		Depth 深度		L	L1	3	H	H1	h	h1		127.30	110.00	85.30	45.00	48.00	43.00	43.90		Bottom 底部[mm]		Position of well center 孔位置[mm]				Theoretical Weight[g]	ID 内径	OD 外径					6	D	P1	P2			4.30	5.60	9.00	9.00		
Plate dimension 板尺寸[mm]		Well dimension 孔尺寸[mm]																																																							
Length 长度		Width 宽度		Height 高度		Depth 深度																																																			
L	L1	3	H	H1	h	h1																																																			
127.30	110.00	85.30	45.00	48.00	43.00	43.90																																																			
Bottom 底部[mm]		Position of well center 孔位置[mm]				Theoretical Weight[g]																																																			
ID 内径	OD 外径																																																								
6	D	P1	P2																																																						
4.30	5.60	9.00	9.00																																																						
Colour	Clear																																																								
Bottle Shape	V type																																																								
Material Characteristics	polypropylene (PP), premium surface																																																								
Product Characteristics	biologically inert																																																								
	no heavy metals according to EC/94/62																																																								
Purity Grade	free of RNase and DNase																																																								
	free of endotoxins																																																								
Shelf - Life	36 months at controlled storage conditions																																																								
Packaging	pcs per bag: 10																																																								
	bag per case: 10																																																								
Storage	Storage Temperature: 2°C-40°C, humidity 40%-80%																																																								
	Keep clean, ventilated and out of direct sunlight																																																								
Approved by PakGent quality department																																																									