


## Technical Datasheet

<b>Article - No:</b>	TDS-CL-MC01	Revision: A0
<b>Remarks</b>	single use only	
Specification		
<b>Product name</b>	CL-MC01, 1layer Cell Factory	
<b>Production conditions</b>	Class 100,000 clean room production	
	No release agent is used	
	Motor molding	
<b>Classification</b>	CE/IVD	
<b>Dimensions</b>		
<b>Cell culture area</b>	632cm <sup>2</sup>	
<b>Material Characteristics</b>	bottle:Polystyrene(PS)	
	cap:High Density Polyethylene(HDPE)	
<b>Product Characteristics</b>	biologically inert	
	no heavy metals according to EC/94/62	
<b>Purity Grade</b>	free of RNase and DNase	
	free of pyrogen	
	free of endotoxins	
	free of ATP	
<b>Shelf - Life</b>	36 months at controlled storage conditions	
<b>Packaging</b>	pcs per case: 6	
<b>Storage</b>	Storage Temperature:2°C-40°C,humidity40%-80%	
	Keep clean,ventilated and out of direct sunlight	
<b>Approved by PakGent quality department</b>		