

Technical Datasheet

Article - No:	TDS-CL-DT035	Revision: A0
Remarks	single use only	
Specification		
Product name	CL-DT035,35mm,Cell Culture Dish,TC treated	
Production conditions	Class 10,000 clean room production	
	No release agent is used	
Classification	CE/IVD	
Dimensions	<p>The drawing shows three views of the dish: a top view (a circle), a side view (a cylinder with a flange), and a bottom view (a circle with a central hole). Dimensions are provided in millimeters: the top diameter is 37.95, the inner diameter of the well is 36.82, the bottom diameter is 34.78, the height of the dish is 12.38, and the diameter of the central hole is 38.79.</p>	
Colour	Clear	
Nominal Capacity	35mm	
Material Characteristics	polystyrene (PS), premium surface	
Sterile	Yes	
Product Characteristics	biologically inert	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogen	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
Shelf - Life	36 months at controlled storage conditions	
	free of ATP	
Packaging	pcs per bag: 20	
	bag per case: 25	
Storage	Storage Temperature:20°C-25°C,humidity40%-60%	
	Keep clean,ventilated and out of direct sunlight	
Approved by PakGent quality department		