


Technical Datasheet

Article - No:	TDS-RT050HF-L	Revision: A0
Remarks	single use only	
Specification		
Product name	RT050HF-L, 50 μ L Hamilton robotic filter tips, low retention	
Production conditions	Class 100,000 clean room production	
	No release agent is used	
	Motor molding	
Classification	CE/IVD	
Dimensions	 <p style="text-align: center;">A - A</p>	
Colour	Clear	
Nominal Capacity	1-50 μ L	
Material Characteristics	polypropylene (PP), premium surface	
	filter: high density polyethylene (PE)	
Product Characteristics	biologically inert	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogen	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
	free of ATP	
Match the model	Star [®]	
	Starlet	
	Nimbus [®]	
Shelf - Life	36 months at controlled storage conditions	
Packaging	pcs per rack: 96	
	rack per case: 30	
Storage	Storage Temperature: 2°C-40°C, humidity 40%-80%	
	Keep clean, ventilated and out of direct sunlight	
Approved by PakGent quality department		