

Technical Datasheet

Article - No:	TDS-SMB125T	Revision: A0
Remarks	single use only	
Specification		
Product name	SMB125T,125ml,Square Media Bottle	
Production conditions	Class 100,000 clean room production	
	No release agent is used	
Classification	CE/IVD	
Dimensions		
Colour	Clear	
Nominal Capacity	125ml	
Material Characteristics	bottle:Polyethylene terephthalate(PET)	
	cap:High Density Polyethylene(HDPE)	
Sterile	Yes	
Product Characteristics	biologically inert	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogen	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
	free of ATP	
Shelf - Life	36 months at controlled storage conditions	
Packaging	pcs per case: 96	
Storage	Storage Temperature:20°C-25°C,humidity40%-60%	
	Keep clean,ventilated and out of direct sunlight	
Approved by PakGent quality department		