


Technical Datasheet

Article - No:	TDS-CL-FT075V	Revision: A0
Remarks	single use only	
Specification		
Product name	CL-FT075V,75cm ² ,Cell Culture Flask,TC Treated,Vent Cap	
Production conditions	Class 10,000 clean room production	
	No release agent is used	
Classification	CE/IVD	
Dimensions		
Colour	Clear	
Graduation	moulded rings	
Nominal Capacity	75cm ²	
Material Characteristics	bottle:polystyrene (PS)	
	cap:High Density Polyethylene(HDPE)	
Sterile	Yes	
Product Characteristics	biologically inert	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogen	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
	free of ATP	
Shelf - Life	36 months at controlled storage conditions	
Packaging	pcs per bag: 5	
	bag per case: 20	
Storage	Storage Temperature:20°C-25°C,humidity40%-60%	
	Keep clean,ventilated and out of direct sunlight	
Approved by PakGent quality department		