

Technical Datasheet

Article - No:	TDS-CL-D150	Revision: A0
Remarks	single use only	
Specification		
Product name	CL-D150,150mm,Cell Culture Dish	
Production conditions	Class 10,000 clean room production	
	No release agent is used	
Classification	CE/IVD	
Dimensions	<p>The drawing shows three views of the dish: a top view, a side view, and a bottom view. The top view is a circle with a diameter of 145.36 mm. The side view shows a depth of 14.9 mm, a rim thickness of 1.1 mm, and a base diameter of 137.55 mm. The bottom view shows a diameter of 145.36 mm and a small circular feature in the center.</p>	
Colour	Clear	
Nominal Capacity	150mm	
Material Characteristics	polystyrene (PS), premium surface	
Sterile	Yes	
Product Characteristics	biologically inert	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogen	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
	free of ATP	
Shelf - Life	36 months at controlled storage conditions	
Packaging	pcs per bag: 5	
	bag per case: 20	
Storage	Storage Temperature:20°C-25°C,humidity40%-60%	
	Keep clean,ventilated and out of direct sunlight	
Approved by PakGent quality department		