

Technical Datasheet

Article - No:	TDS-RT200TF-L	Revision: A0
Remarks	single use only	
Specification		
Product name	RT200TF-L,200µL Tecan robotic filter tips,low retention	
Production conditions	Class 100,000 clean room production	
	No release agent is used	
	Motor molding	
Classification	CE/IVD	
Dimensions		
Colour	Clear	
Nominal Capacity	1-200µL	
Material Characteristics	polypropylene (PP),premium surface	
	filter:high density polyethylene (PE)	
Product Characteristics	biologically inert	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogen	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
Match the model	free of ATP	
	Genesis Freedom®	
	Freedom EVO®	
Shelf - Life	Miniprep with LiHa®	
	36 months at controlled storage conditions	
Packaging	pcs per rack: 96	
	rack per case: 30	
Storage	Storage Temperature:2°C-40°C,humidity40%-80%	
	Keep clean,ventilated and out of direct sunlight	
Approved by PakGent quality department		