

Technical Datasheet

Article - No:	TDS-RT20R-L	Revision: A0
Remarks	single use only	
Specification		
Product name	RT20R-L,20µL Rainin LTS Pipette tips,low retention	
Production conditions	Class 100,000 clean room production	
	No release agent is used	
	Motor molding	
Classification	CE/IVD	
Dimensions		
Colour	Clear	
Graduation	moulded rings	
Nominal Capacity	0.1-20µL	
Material Characteristics	polypropylene (PP), premium surface	
Product Characteristics	biologically inert	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogen	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
	free of ATP	
Shelf - Life	36 months at controlled storage conditions	
Packaging	pcs per box: 96	
	box per case: 50	
Storage	Storage Temperature:2°C-40°C,humidity40%-80%	
	Keep clean,ventilated and out of direct sunlight	
Approved by PakGent quality department		