




Pipette tips



Catalog

+86 (0) 512 6670 0017
info@pakgentbio.com
www.pakgentbio.com



Pipette tips



All of our tips are tested by an independent laboratory and certified to be RNase and DNase free. Sterile products are certified to be free from RNase, DNase, DNA, pyrogens, bioburden, PCR inhibitors, and verified sterile.

Fit to Brand Pipettes


[Eppendorf](#) [Thermo](#) [Biohit](#) [Glison](#) [Brand](#) [Rainin](#) [DLAB](#) [IKA.....](#)

Low Retention Pipette tips

- It is more suitable for taking liquids with high viscosity, high oil property and detergent.
- It is recommended to use a low retention tip for precious liquids.



Non Low Retention Low Retention




Larger Volume Pipette Tips

- Larger volume, but more accurate.
- It can reduce or eliminate the use of fragile glass and plastic serological pipettes. Hands-free tip loading is now possible when using 5 ml and 10 ml pipettes.

Filter Tips


- The filter tips is especially suitable for PCR and experiments with aerosol pollution.
- Commonly used in cell culture, Bacteriology, Virology, molecular biology and liquid transfer of radioactive samples.



Refill Pack Pipette Tips


- Patented R&D blister box packaging and product pre-sterilization ensure convenient use and space saving, The unused pipette tips shall be kept in a relatively closed space to ensure the sterilization effect.

- The number of stiffeners is increased, and the angle is 45 degrees, which increases the material strength and better protects the product. Suitable for mass packaging and various transportation environments.



LTS Tips (Light Touch Pipette Tips)

- LTS pipette tips users can achieve reliable and consistent fitting effect without excessive force. LTS pipette tips has well-designed small sealing area and front stop point design. The front stop design allows you to know the exact position of the fitting.



Tip Ejection Forces

	0	1	2	3	4kg
With LTS	0.6kg				
Without LTS					4kg
Recommended 30% Female Max (Strength Capacity)			2.1kg		
Recommended 30% Male Max (Strength Capacity)				3kg	

Universal Pipette Tips

10uL standard, 32mm



Part NO.	Volume	Filter	Steril	Pack
UFPT-10R	0.1-10 ul	x	✓	96/rack, 50racks/case
UFPT-F-10	0.1-10 ul	✓	✓	96/rack, 50racks/case
UFPT-REF-10	0.1-10 ul	x	✓	96/rack, 10racks/unit, 10units/case
UFPT-10B	0.1-10 ul	x	x	1,000/bag, 40bags/case
UFPT-F-10B	0.1-10 ul	✓	x	1,000/bag, 40bags/case

10uL extra long design, 46mm



Part NO.	Volume	Filter	Steril	Pack
UFPT-10ER	0.1-10 ul	x	✓	96/rack, 50racks/case
UFPT-F-10E	0.1-10 ul	✓	✓	96/rack, 50racks/case
UFPT-10EB	0.1-10 ul	x	x	1,000/bag, 40bags/case
UFPT-F-10EB	0.1-10 ul	✓	x	1,000/bag, 40bags/case

200uL standard wide, 51mm



Part NO.	Volume	Filter	Steril	Pack
UFPT-F-50	1-50 ul	✓	✓	96/rack, 50racks/case
UFPT-F-100	1-100 ul	✓	✓	96/rack, 50racks/case
UFPT-F-200	1-200 ul	✓	✓	96/rack, 50racks/case
UFPT-F-50B	1-50 ul	✓	x	1,000/bag, 20bags/case
UFPT-F-100B	1-100 ul	✓	x	1,000/bag, 20bags/case
UFPT-F-200B	1-200 ul	✓	x	1,000/bag, 20bags/case

1000uL standard, 78mm



Part NO.	Volume	Filter	Steril	Pack
UFPT-1000R	100-1,000 ul	x	✓	96/rack, 50racks/case
UFPT-F-1000	100-1,000 ul	✓	✓	96/rack, 50racks/case
UFPT-1000B	100-1,000 ul	x	x	1,000/bag, 10bags/case
UFPT-F-1000B	100-1,000 ul	✓	x	1,000/bag, 10bags/case

10uL extra long, thin and sharp design, 46mm



Part NO.	Volume	Filter	Steril	Pack
UFPT-10ETR	0.1-10 ul	x	✓	96/rack, 50racks/case
UFPT-F-10ET	0.1-10 ul	✓	✓	96/rack, 50racks/case
UFPT-REF-10ET	0.1-10 ul	x	✓	96/rack, 9racks/unit, 10units/case
UFPT-10ETB	0.1-10 ul	x	x	1,000/bag, 40bags/case
UFPT-F-10ETB	0.1-10 ul	✓	x	1,000/bag, 40bags/case

200uL standard, yellow, 51mm



Part NO.	Volume	Filter	Steril	Pack
UFPT-200R	1-200 ul	x	✓	96/rack, 50racks/case
UFPT-F-20	1-20 ul	✓	✓	96/rack, 50racks/case
UFPT-REF-200	1-200 ul	x	✓	96/rack, 12racks/unit, 10units/case
UFPT-200B	1-200 ul	x	x	1,000/bag, 20bags/case
UFPT-F-20B	1-20 ul	✓	x	1,000/bag, 20bags/case

200uL extra long or 300uL, 60mm



Part NO.	Volume	Filter	Steril	Pack
UFPT-200ER	1-300 ul	x	✓	96/rack, 50racks/case
UFPT-F-100E	1-100 ul	✓	✓	96/rack, 50racks/case
UFPT-F-200E	1-300 ul	✓	✓	96/rack, 50racks/case
UFPT-REF-200E	1-300 ul	x	✓	96/rack, 10racks/unit, 10units/case
UFPT-200EB	1-300 ul	x	x	1,000/bag, 20bags/case
UFPT-F-200EB	1-300 ul	✓	x	1,000/bag, 20bags/case

1000uL extra long or 1250uL, 102mm



Part NO.	Volume	Filter	Steril	Pack
UFPT-1000ER	100-1,000 ul	x	✓	96/rack, 50racks/case
UFPT-F-1000E	100-1,000 ul	✓	✓	96/rack, 50racks/case
UFPT-REF-1000E	100-1,000 ul	x	✓	96/rack, 8racks/unit, 10units/case
UFPT-1000EB	100-1,000 ul	x	x	1,000/bag, 10bags/case
UFPT-F-1000EB	100-1,000 ul	✓	x	1,000/bag, 10bags/case

Super long 200uL, 90mm



Part NO.	Volume	Filter	Steril	Pack
UFPT-200SR	1-200 ul	x	✓	96/rack, 50racks/case
UFPT-F-200S	1-200 ul	✓	✓	96/rack, 50racks/case
UFPT-200SB	1-200 ul	x	x	1,000/bag, 10bags/case
UFPT-F-200SB	1-200 ul	✓	x	1,000/bag, 10bags/case

Low Retention Pipette tips

10uL standard, 32mm



Part NO.	Volume	Filter	Steril	Pack
UFPT-10R-L	0.1-10 ul	x	✓	96/rack, 50racks/case
UFPT-F-10L	0.1-10 ul	✓	✓	96/rack, 50racks/case
UFPT-REF-10L	0.1-10 ul	x	✓	96/rack, 10racks/unit, 10units/case
UFPT-10B-L	0.1-10 ul	x	x	1,000/bag, 40bags/case
UFPT-F-10B-L	0.1-10 ul	✓	x	1,000/bag, 40bags/case

10uL extra long design, 46mm



Part NO.	Volume	Filter	Steril	Pack
UFPT-10ER-L	0.1-10 ul	x	✓	96/rack, 50racks/case
UFPT-F-10E-L	0.1-10 ul	✓	✓	96/rack, 50racks/case
UFPT-10EB-L	0.1-10 ul	x	x	1,000/bag, 40bags/case
UFPT-F-10EB-L	0.1-10 ul	✓	x	1,000/bag, 40bags/case

200uL standard wide, 51mm



Part NO.	Volume	Filter	Steril	Pack
UFPT-F-50L	1-50 ul	✓	✓	96/rack, 50racks/case
UFPT-F-100L	1-100 ul	✓	✓	96/rack, 50racks/case
UFPT-F-200L	1-200 ul	✓	✓	96/rack, 50racks/case
UFPT-F-50B-L	1-50 ul	✓	x	1,000/bag, 20bags/case
UFPT-F-100B-L	1-100 ul	✓	x	1,000/bag, 20bags/case
UFPT-F-200B-L	1-200 ul	✓	x	1,000/bag, 20bags/case

Medium long 1000uL, 88mm



Part NO.	Volume	Filter	Steril	Pack
UFPT-1000MR	100-1,000 ul	x	✓	96/rack, 50racks/case
UFPT-F-1000M	100-1,000 ul	✓	✓	96/rack, 50racks/case
UFPT-1000MB	100-1,000 ul	x	x	1,000/bag, 10bags/case
UFPT-F-1000MB	100-1,000 ul	✓	x	1,000/bag, 10bags/case

10uL extra long, thin and sharp design, 46mm



Part NO.	Volume	Filter	Steril	Pack
UFPT-10ETR-L	0.1-10 ul	x	✓	96/rack, 50racks/case
UFPT-F-10ET-L	0.1-10 ul	✓	✓	96/rack, 50racks/case
UFPT-REF-10ET-L	0.1-10 ul	x	✓	96/rack, 9racks/unit, 10units/case
UFPT-10ETB-L	0.1-10 ul	x	x	1,000/bag, 40bags/case
UFPT-F-10ETB-L	0.1-10 ul	✓	x	1,000/bag, 40bags/case

200uL standard, yellow, 51mm



Part NO.	Volume	Filter	Steril	Pack
UFPT-200R-L	1-200 ul	x	✓	96/rack, 50racks/case
UFPT-F-20L	1-20 ul	✓	✓	96/rack, 50racks/case
UFPT-REF-200L	1-200 ul	x	✓	96/rack, 12racks/unit, 10units/case
UFPT-200B-L	1-200 ul	x	x	1,000/bag, 20bags/case
UFPT-F-20B-L	1-20 ul	✓	x	1,000/bag, 20bags/case

200uL extra long or 300uL, 60mm



Part NO.	Volume	Filter	Steril	Pack
UFPT-200ER-L	1-300 ul	x	✓	96/rack, 50racks/case
UFPT-F-100E-L	1-100 ul	✓	✓	96/rack, 50racks/case
UFPT-F-200E-L	1-300 ul	✓	✓	96/rack, 50racks/case
UFPT-REF-200E-L	1-300 ul	x	✓	96/rack, 10racks/unit, 10units/case
UFPT-200EB-L	1-300 ul	x	x	1,000/bag, 20bags/case
UFPT-F-200EB-L	1-300 ul	✓	x	1,000/bag, 20bags/case

1000uL standard, 78mm



Part NO.	Volume	Filter	Steril	Pack
UFPT-1000R-L	100-1,000 ul	x	✓	96/rack, 50racks/case
UFPT-F-1000-L	100-1,000 ul	✓	✓	96/rack, 50racks/case
UFPT-1000B-L	100-1,000 ul	x	x	1,000/bag, 10bags/case
UFPT-F-1000B-L	100-1,000 ul	✓	x	1,000/bag, 10bags/case

Super long 200uL, 90mm



Part NO.	Volume	Filter	Steril	Pack
UFPT-200SR-L	1-200 ul	x	✓	96/rack, 50racks/case
UFPT-F-200S-L	1-200 ul	✓	✓	96/rack, 50racks/case
UFPT-200SB-L	1-200 ul	x	x	1,000/bag, 10bags/case
UFPT-F-200SB-L	1-200 ul	✓	x	1,000/bag, 10bags/case

Larger Volume Pipette Tips

5ml, clear, fit Eppendorf, Biohit, DLAB Pipette



Part NO.	Filter	Steril	Pack
PT-5000R-E	x	✓	24/rack, 30racks/case
PT-5000RF-E	✓	✓	24/rack, 30racks/case
PT-5000B-E	x	x	100/bag, 20bags/case
PT-5000BF-E	✓	x	100/bag, 20bags/case

10ml, clear, fit Eppendorf, Biohit, DLAB Pipette



Part NO.	Filter	Steril	Pack
PT-10000R-E	x	✓	24/rack, 40racks/case
PT-10000RF-E	✓	✓	24/rack, 40racks/case
PT-10000B-E	x	x	100/bag, 10bags/case
PT-10000BF-E	✓	x	100/bag, 10bags/case

1000uL extra long or 1250uL, 102mm



Part NO.	Volume	Filter	Steril	Pack
UFPT-1000ER-L	100-1,000 ul	x	✓	96/rack, 50racks/case
UFPT-F-1000E-L	100-1,000 ul	✓	✓	96/rack, 50racks/case
UFPT-REF-1000E-L	100-1,000 ul	x	✓	96/rack, 8racks/unit, 10units/case
UFPT-1000EB-L	100-1,000 ul	x	x	1,000/bag, 10bags/case
UFPT-F-1000EB-L	100-1,000 ul	✓	x	1,000/bag, 10bags/case

Medium long 1000uL, 88mm



Part NO.	Volume	Filter	Steril	Pack
UFPT-1000MR-L	100-1,000 ul	x	✓	96/rack, 50racks/case
UFPT-F-1000M-L	100-1,000 ul	✓	✓	96/rack, 50racks/case
UFPT-1000MB-L	100-1,000 ul	x	x	1,000/bag, 10bags/case
UFPT-F-1000MB-L	100-1,000 ul	✓	x	1,000/bag, 10bags/case

5ml, clear, fit Thermo, Brand Pipette



Part NO.	Filter	Steril	Pack
PT-5000R-T	x	✓	40/rack, 40racks/case
PT-5000RF-T	✓	✓	40/rack, 40racks/case
PT-5000B-T	x	x	100/bag, 20bags/case
PT-5000BF-T	✓	x	100/bag, 20bags/case

10ml, clear, fit Thermo, Brand Pipette



Part NO.	Filter	Steril	Pack
PT-10000R-T	x	✓	24/rack, 40racks/case
PT-10000RF-T	✓	✓	24/rack, 40racks/case
PT-10000B-T	x	x	100/bag, 10bags/case
PT-10000BF-T	✓	x	100/bag, 10bags/case

LTS Tips (Light Touch Pipette Tips)

20ul Fit LTS Pipette, sterile



Part NO.	Filter	Low Retention	Steril	Pack
RT20R	x	x	✓	96/rack, 50racks/case
RT20R-L	x	✓	✓	96/rack, 50racks/case
RT20RF	✓	x	✓	96/rack, 50racks/case
RT20RF-L	✓	✓	✓	96/rack, 50racks/case
RT20B	x	x	x	1,000/bag, 40bags/case
RT20B-L	x	✓	x	1,000/bag, 40bags/case
RT20-REF	x	x	✓	96/rack, 10racks/unit, 10units/case
RT20-REF-L	x	✓	✓	96/rack, 10racks/unit, 10units/case

300ul Fit LTS Pipette, sterile



Part NO.	Filter	Low Retention	Steril	Pack
RT300R	x	x	✓	96/rack, 50racks/case
RT300R-L	x	✓	✓	96/rack, 50racks/case
RT300RF	✓	x	✓	96/rack, 50racks/case
RT300RF-L	✓	✓	✓	96/rack, 50racks/case
RT300B	x	x	x	1,000/bag, 20bags/case
RT300B-L	x	✓	x	1,000/bag, 20bags/case
RT300-REF	x	x	✓	96/rack, 8racks/unit, 10units/case
RT300-REF-L	x	✓	✓	96/rack, 8racks/unit, 10units/case

1200ul Fit LTS Pipette, sterile



Part NO.	Filter	Low Retention	Steril	Pack
RT1200R	x	x	✓	96/rack, 50racks/case
RT1200R-L	x	✓	✓	96/rack, 50racks/case
RT1200RF	✓	x	✓	96/rack, 50racks/case
RT1200RF-L	✓	✓	✓	96/rack, 50racks/case
RT1200B	x	x	x	1,000/bag, 10bags/case
RT1200B-L	x	✓	x	1,000/bag, 10bags/case

200ul Fit LTS Pipette, sterile




Part NO.	Filter	Low Retention	Steril	Pack
RT200R	x	x	✓	96/rack, 50racks/case
RT200R-L	x	✓	✓	96/rack, 50racks/case
RT200RF	✓	x	✓	96/rack, 50racks/case
RT200RF-L	✓	✓	✓	96/rack, 50racks/case
RT200B	x	x	x	1,000/bag, 20bags/case
RT200B-L	x	✓	x	1,000/bag, 20bags/case
RT200-REF	x	x	✓	96/rack, 10racks/unit, 10units/case
RT200-REF-L	x	✓	✓	96/rack, 10racks/unit, 10units/case

1000ul Fit LTS Pipette, sterile



Part NO.	Filter	Low Retention	Steril	Pack
RT1000R	x	x	✓	96/rack, 50racks/case
RT1000R-L	x	✓	✓	96/rack, 50racks/case
RT1000RF	✓	x	✓	96/rack, 50racks/case
RT1000RF-L	✓	✓	✓	96/rack, 50racks/case
RT1000B	x	x	x	1,000/bag, 10bags/case
RT1000B-L	x	✓	x	1,000/bag, 10bags/case
RT1000-REF	x	x	✓	96/rack, 8racks/unit, 10units/case
RT1000-REF-L	x	✓	✓	96/rack, 8racks/unit, 10units/case



Rack / Bulk / Refill Pack Available