

Technical Datasheet

Article - No:	TDS-UFPT-REF-200-L	Revision: A0
Remarks	single use only	
Specification		
Product name	UFPT-REF-200-L,200µL Pipette tips,Refull pack,low retention	
Production conditions	Class 100,000 clean room production	
	No release agent is used	
	Motor molding	
Classification	CE/IVD	
Dimensions		
Colour	Clear	
Graduation	moulded rings	
Nominal Capacity	1-200µL	
Sterile	Yes	
Material Characteristics	polypropylene(PP),premium surface	
Product Characteristics	pack:Polyethylene terephthalate(PET)	
	biologically inert	
Purity Grade	no heavy metals according to EC/94/62	
	free of RNase and DNase	
	free of pyrogen	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
Shelf - Life	free of ATP	
	36 months at controlled storage conditions	
Packaging	pcs per rack: 96	
	rack per pack: 12	
	pack per case: 10	
Storage	Storage Temperature:2°C-40°C,humidity40%-80%	
	Keep clean,ventilated and out of direct sunlight	
Approved by PakGent quality department		