


Technical Datasheet

Article - No:	TDS-EFG-500V	Revision: A0
Remarks	single use only	
Specification		
Product name	EFG-500V, 500ml, Erlenmeyer Flask, Vent Cap, Flat bottom	
Production conditions	Class 100,000 clean room production	
	No release agent is used	
Classification	CE/IVD	
Dimensions		
Colour	Clear	
Graduation	moulded rings	
Nominal Capacity	500ml	
Material Characteristics	bottle: Polyethylene Terephthalate Glycol (PETG)	
	cap: High Density Polyethylene (HDPE)	
Sterile	Yes	
Product Characteristics	biologically inert	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogen	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
	free of ATP	
Shelf - Life	36 months at controlled storage conditions	
Packaging	pcs per case: 12	
Storage	Storage Temperature: 20°C-25°C, humidity 40%-60%	
	Keep clean, ventilated and out of direct sunlight	
Approved by PakGent quality department		