

Technical Datasheet

Article - No:	TDS-SPTS50	Revision: A0
Remarks	single use only	
Specification		
Product name	SPTS50,50ML Serological Pipette,Individual packaging	
Production conditions	Class 100,000 clean room production	
	No release agent is used	
	Welding moulding	
Classification	CE/IVD	
Dimensions	<p style="text-align: center;"> 354 ± 2 40 ± 1 26 ± 1 18 ± 0.1 16.07 ± 0.1 32 ± 0.5 5.5 ± 0.5 </p>	
Colour	Clear	
Graduation	Printing	
Nominal Capacity	50ML	
Sterile	Yes	
Material Characteristics	Polystyrene (PS), premium surface	
	filter:high density polyethylene (PE)	
Product Characteristics	biologically inert	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogen	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
	free of ATP	
Shelf - Life	36 months at controlled storage conditions	
Packaging	pcs per bag: 25	
	bag per case: 8	
Storage	Storage Temperature:2°C-40°C,humidity40%-80%	
	Keep clean,ventilated and out of direct sunlight	
Approved by PakGent quality department		