

Technical Datasheet

Article - No:	TDS-CL-DT060	Revision: A0
Remarks	single use only	
	Specification	Specification
Product name	CL-DT060,60mm,Cell Culture Dish,TC treated	
Production conditions	Class 10,000 clean room production	
	No release agent is used	
Classification	CE/IVD	
Dimensions		
Colour	Clear	
Nominal Capacity	60mm	
Material Characteristics	polystyrene (PS), premium surface	
Sterile	Yes	
Product Characteristics	biologically inert	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogen	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
Shelf - Life	36 months at controlled storage conditions	
	free of ATP	
Packaging	pcs per bag: 10	
	bag per case: 50	
Storage	Storage Temperature:20°C-25°C,humidity40%-60%	
	Keep clean,ventilated and out of direct sunlight	
Approved by PakGent quality department		