

Technical Datasheet

| | | |
|---|--|--------------|
| Article - No: | TDS-RT1200B | Revision: A0 |
| Remarks | single use only | |
| Specification | | |
| Product name | RT1200B, 1200µL Rainin LTS Pipette tips | |
| Production conditions | Class 100,000 clean room production | |
| | No release agent is used | |
| | Motor molding | |
| Classification | CE/IVD | |
| Dimensions | <p>The drawing shows a pipette tip with a total length of 86.6. The base has a diameter of 8.1. The neck has a diameter of 8.8. The tip is labeled 20ul and has a diameter of 1.34±0.05. The nozzle has a diameter of 0.95±0.05.</p> | |
| Colour | Clear | |
| Graduation | moulded rings | |
| Nominal Capacity | 100-1200µL | |
| Material Characteristics | polypropylene (PP), premium surface | |
| Product Characteristics | biologically inert | |
| | no heavy metals according to EC/94/62 | |
| Purity Grade | free of RNase and DNase | |
| | free of pyrogen | |
| | free of PCR inhibitors | |
| | free of human DNA | |
| | free of endotoxins | |
| Shelf - Life | 36 months at controlled storage conditions | |
| Packaging | pcs per bag: 1000 | |
| | bag per case: 10 | |
| Storage | Storage Temperature: 2°C-40°C, humidity 40%-80% | |
| | Keep clean, ventilated and out of direct sunlight | |
| Approved by PakGent quality department | | |