

## Technical Datasheet

<b>Article - No:</b>	TDS-CL-F225V	Revision: A0
<b>Remarks</b>	single use only	
<b>Specification</b>		
<b>Product name</b>	CL-F225V, 225cm <sup>2</sup> , Cell Culture Flask, Vent cap	
<b>Production conditions</b>	Class 10,000 clean room production	
	No release agent is used	
<b>Classification</b>	CE/IVD	
<b>Dimensions</b>		
<b>Colour</b>	Clear	
<b>Graduation</b>	moulded rings	
<b>Nominal Capacity</b>	225cm <sup>2</sup>	
<b>Material Characteristics</b>	bottle: polystyrene (PS)	
	cap: High Density Polyethylene (HDPE)	
<b>Sterile</b>	Yes	
<b>Product Characteristics</b>	biologically inert	
	physiologically harmless	
	no heavy metals according to EC/94/62	
<b>Purity Grade</b>	free of RNase and DNase	
	free of pyrogen	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
	free of ATP	
<b>Shelf - Life</b>	36 months at controlled storage conditions	
<b>Packaging</b>	pcs per bag: 5	
	bag per case: 5	
<b>Storage</b>	Storage Temperature: 20°C-25°C, humidity 40%-60%	
	Keep clean, ventilated and out of direct sunlight	
<b>Approved by PakGent quality department</b>		