

Technical Datasheet

Article - No:	TDS-NMPB008	Revision: A0
Remarks	single use only	
Specification		
Product name	NMPB008,8ml Packaging Bottle,HDPE	
Production conditions	Class 100,000 clean room production	
	No release agent is used	
Classification	CE/IVDD	
Dimensions		
Colour	Clear	
Nominal Capacity	8ml	
Material Characteristics	bottle:High Density Polyethylene(HDPE)	
	cap:polypropylene (PP)	
Performance parameters	Negative pressure:-70kpa~-90kpa	
	Ethanol loss rate: $\leq 0.5\%$	
	water loss rate: $\leq 0.2\%$	
Purity Grade	free of RNase and DNase	
	free of pyrogen	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
Shelf - Life	36 months at controlled storage conditions	
Packaging	pcs per case:2000	
Storage	dry condition:40%-60% rel. humidity	
Approved by PakGent quality department		