


## Technical Datasheet

|   |  |              |
|---|--|--------------|
| <b>Article - No:</b>                          | TDS-PT-5000BF-T  | Revision: A0 |
| <b>Remarks</b>                                | single use only  |              |
| Specification                                 |  |              |
| <b>Product name</b>                           | PT-5000BF-T,5mL filter tips  |              |
| <b>Production conditions</b>                  | Class 100,000 clean room production  |              |
|   | No release agent is used   |              |
|   | Motor molding  |              |
| <b>Classification</b>                         | CE/IVD   |              |
| <b>Dimensions</b>                             |  |              |
| <b>Colour</b>                                 | Clear  |              |
| <b>Graduation</b>                             | moulded rings  |              |
| <b>Nominal Capacity</b>                       | 5mL  |              |
| <b>Material Characteristics</b>               | polypropylene(PP),premium surface  |              |
|   | filter:high density polyethylene(PE)   |              |
| <b>Product Characteristics</b>                | biologically inert   |              |
|   | no heavy metals according to EC/94/62  |              |
| <b>Purity Grade</b>                           | free of RNase and DNase  |              |
|   | free of pyrogen  |              |
|   | free of PCR inhibitors   |              |
|   | free of human DNA  |              |
| <b>Match the model</b>                        | Thermo   |              |
|   | Brand  |              |
| <b>Shelf - Life</b>                           | 36 months at controlled storage conditions   |              |
| <b>Packaging</b>                              | pcs per bag: 100   |              |
|   | bag per case: 20   |              |
| <b>Storage</b>                                | Storage Temperature:2°C-40°C,humidity40%-80%                                       |              |
|   | Keep clean,ventilated and out of direct sunlight                                   |              |
| <b>Approved by PakGent quality department</b> |  |              |