

Technical Datasheet

| | | |
|---|--|--------------|
| Article - No: | TDS-STSS-509AM | Revision: A0 |
| Remarks | single use only | |
| Specification | | |
| Product name | STSS-509AM,0.5ml Screw Cap Tubes | |
| Production conditions | Class 100,000 clean room production | |
| | No release agent is used | |
| Classification | CE/IVDD | |
| Dimensions | | |
| Colour | Amber tube,Amber cap | |
| Nominal Capacity | 0.5ml | |
| Material Characteristics | polypropylene (PP) | |
| Performance parameters | Storage Temperature: -80°C~121°C | |
| Product Characteristics | biologically inert | |
| | no heavy metals according to EC/94/62 | |
| Purity Grade | free of RNase and DNase | |
| | free of pyrogen | |
| | free of PCR inhibitors | |
| | free of human DNA | |
| | free of endotoxins | |
| Shelf - Life | 36 months at controlled storage conditions | |
| Packaging | pcs per bag:500 | |
| | bag per case:10 | |
| Storage | dry condition:40%-60% rel. humidity | |
| Approved by PakGent quality department | | |