


Technical Datasheet

Article - No:	TDS-PT-10000BF-E	Revision: A0
Remarks	single use only	
Specification		
Product name	PT-10000BF-E, 10mL filter tips	
Production conditions	Class 100,000 clean room production	
	No release agent is used	
	Motor molding	
Classification	CE/IVD	
Dimensions		
Colour	Clear	
Graduation	moulded rings	
Nominal Capacity	10mL	
Material Characteristics	polypropylene(PP), premium surface	
	filter: high density polyethylene(PE)	
Product Characteristics	biologically inert	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogen	
	free of PCR inhibitors	
	free of human DNA	
Match the model	free of endotoxins	
	Eppendorf	
	Biohit	
Shelf - Life	DLAB	
	36 months at controlled storage conditions	
Packaging	pcs per bag: 100	
	bag per case: 10	
Storage	Storage Temperature: 2°C-40°C, humidity 40%-80%	
	Keep clean, ventilated and out of direct sunlight	
Approved by PakGent quality department		