

Technical Datasheet

Article - No:	TDS-SPTS10A	Revision: A0
Remarks	single use only	
Specification		
Product name	SPTS10A,10ML Serological Pipette,Individual packaging	
Production conditions	Class 100,000 clean room production	
	No release agent is used	
	Welding moulding	
Classification	CE/IVD	
Dimensions		
Colour	Clear	
Graduation	Printing	
Nominal Capacity	10ML	
Sterile	Yes	
Material Characteristics	Polystyrene (PS), premium surface	
	filter:high density polyethylene (PE)	
Product Characteristics	biologically inert	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogen	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
Shelf - Life	free of ATP	
	36 months at controlled storage conditions	
Packaging	pcs per bag: 50	
	bag per case: 10	
Storage	Storage Temperature:2°C-40°C,humidity40%-80%	
	Keep clean,ventilated and out of direct sunlight	
Approved by PakGent quality department		