


Technical Datasheet

Article - No:	TDS-NMPB060PA	Revision: A0																						
Remarks	single use only																							
Specification																								
Product name	NMPB060PA,60ml Packaging Bottle,PP,Amber																							
Production conditions	Class 100,000 clean room production																							
	No release agent is used																							
Classification	CE/IVDD																							
Dimensions																								
	<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th>location name</th> <th>material</th> <th>O.D. D1 (mm)</th> <th>Height H1 (mm)</th> <th>I.D. D2 (mm)</th> <th>Cap O.D. D3 (mm)</th> <th>Cap H H2 (mm)</th> <th>Total H (mm)</th> </tr> </thead> <tbody> <tr> <td rowspan="2" style="text-align: center;">60ml</td> <td style="text-align: center;">PP</td> <td style="text-align: center;">φ39.0</td> <td style="text-align: center;">84.3</td> <td style="text-align: center;">φ13.8</td> <td style="text-align: center;">φ24.2</td> <td style="text-align: center;">14.5</td> <td style="text-align: center;">85.5</td> </tr> <tr> <td style="text-align: center;">PE</td> <td style="text-align: center;">φ39.0</td> <td style="text-align: center;">84.3</td> <td style="text-align: center;">φ13.8</td> <td style="text-align: center;">φ24.2</td> <td style="text-align: center;">14.5</td> <td style="text-align: center;">85.5</td> </tr> </tbody> </table>		location name	material	O.D. D1 (mm)	Height H1 (mm)	I.D. D2 (mm)	Cap O.D. D3 (mm)	Cap H H2 (mm)	Total H (mm)	60ml	PP	φ39.0	84.3	φ13.8	φ24.2	14.5	85.5	PE	φ39.0	84.3	φ13.8	φ24.2	14.5
location name	material	O.D. D1 (mm)	Height H1 (mm)	I.D. D2 (mm)	Cap O.D. D3 (mm)	Cap H H2 (mm)	Total H (mm)																	
60ml	PP	φ39.0	84.3	φ13.8	φ24.2	14.5	85.5																	
	PE	φ39.0	84.3	φ13.8	φ24.2	14.5	85.5																	
Colour	Amber																							
Nominal Capacity	60ml																							
Material Characteristics	bottle:polypropylene (PP)																							
	cap:High Density Polyethylene(HDPE)																							
Performance parameters	Negative pressure:-70kpa~-95kpa																							
	Ethanol loss rate: ≤0.5%																							
	water loss rate: ≤0.2%																							
Purity Grade	free of RNase and DNase																							
	free of pyrogen																							
	free of PCR inhibitors																							
	free of human DNA																							
	free of endotoxins																							
Shelf - Life	36 months at controlled storage conditions																							
Packaging	pcs per case:800																							
Storage	dry condition:40%-60% rel. humidity																							
Approved by PakGent quality department																								