

## Technical Datasheet

|   |   |              |
|---|---|--------------|
| <b>Article - No:</b>                          | TDS-PT-10000R-T                                   | Revision: A0 |
| <b>Remarks</b>                                | single use only                                   |              |
| Specification                                 |   |              |
| <b>Product name</b>                           | PT-10000R-T, 10mL Pipette tips                    |              |
| <b>Production conditions</b>                  | Class 100,000 clean room production               |              |
|   | No release agent is used                          |              |
|   | Motor molding                                     |              |
| <b>Classification</b>                         | CE/IVD  |              |
| <b>Dimensions</b>                             |   |              |
| <b>Colour</b>                                 | Clear   |              |
| <b>Graduation</b>                             | moulded rings                                     |              |
| <b>Nominal Capacity</b>                       | 10mL  |              |
| <b>Material Characteristics</b>               | polypropylene(PP), premium surface                |              |
| <b>Product Characteristics</b>                | biologically inert                                |              |
|   | no heavy metals according to EC/94/62             |              |
| <b>Purity Grade</b>                           | free of RNase and DNase                           |              |
|   | free of pyrogen                                   |              |
|   | free of PCR inhibitors                            |              |
|   | free of human DNA                                 |              |
|   | free of endotoxins                                |              |
| <b>Match the model</b>                        | Thermo  |              |
|   | Brand   |              |
| <b>Shelf - Life</b>                           | 36 months at controlled storage conditions        |              |
| <b>Packaging</b>                              | pcs per rack: 24                                  |              |
|   | rack per case: 40                                 |              |
| <b>Storage</b>                                | Storage Temperature: 2°C-40°C, humidity 40%-80%   |              |
|   | Keep clean, ventilated and out of direct sunlight |              |
| <b>Approved by PakGent quality department</b> |   |              |