


## Technical Datasheet

<b>Article - No:</b>	TDS-UFPT-REF-1000E	Revision: A0
<b>Remarks</b>	single use only	
<b>Specification</b>		
<b>Product name</b>	UFPT-REF-1000E, 1000µL Extra long Pipette tips, Refull pack	
<b>Production conditions</b>	Class 100,000 clean room production	
	No release agent is used	
	Motor molding	
<b>Classification</b>	CE/IVD	
<b>Dimensions</b>	 <p>The drawing shows a side view of a pipette tip. The total length is 102<sup>+0</sup><sub>-1</sub> mm. The diameter at the base is 8.9 mm, and at the top is 7.9 mm. The tip has a conical shape with a diameter of 1.3±0.05 mm at the very tip and 0.81±0.05 mm at the base of the tip. Graduations are marked at 1000µL, 500µL, and 200µL.</p>	
<b>Colour</b>	Blue	
<b>Graduation</b>	moulded rings	
<b>Nominal Capacity</b>	100-1200µL	
<b>Sterile</b>	Yes	
<b>Material Characteristics</b>	polypropylene(PP), premium surface	
	pack: Polyethylene terephthalate(PET)	
<b>Product Characteristics</b>	biologically inert	
	no heavy metals according to EC/94/62	
	free of RNase and DNase	
<b>Purity Grade</b>	free of pyrogen	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
	free of ATP	
<b>Shelf - Life</b>	36 months at controlled storage conditions	
<b>Packaging</b>	pcs per rack: 96	
	rack per pack: 8	
	pack per case: 10	
<b>Storage</b>	Storage Temperature: 2°C-40°C, humidity 40%-80%	
	Keep clean, ventilated and out of direct sunlight	
<b>Approved by PakGent quality department</b>		