


## Technical Datasheet

<b>Article - No:</b>	TDS-UFPT-F-20B	Revision: A0
<b>Remarks</b>	single use only	
Specification		
<b>Product name</b>	UFPT-F-20B,20µL filter tips	
<b>Production conditions</b>	Class 100,000 clean room production	
	No release agent is used	
	Motor molding	
<b>Classification</b>	CE/IVD	
<b>Dimensions</b>		
<b>Colour</b>	Yellow	
<b>Graduation</b>	moulded rings	
<b>Nominal Capacity</b>	1-20µL	
<b>Material Characteristics</b>	polypropylene(PP),premium surface	
	filter:high density polyethylene(PE)	
<b>Product Characteristics</b>	biologically inert	
	no heavy metals according to EC/94/62	
<b>Purity Grade</b>	free of RNase and DNase	
	free of pyrogen	
	free of PCR inhibitors	
	free of human DNA	
<b>Shelf - Life</b>	free of endotoxins	
	36 months at controlled storage conditions	
<b>Packaging</b>	pcs per bag: 1000	
	bag per case: 20	
<b>Storage</b>	Storage Temperature:2°C-40°C,humidity40%-80%	
	Keep clean, ventilated and out of direct sunlight	
<b>Approved by PakGent quality department</b>		