


Technical Datasheet

Article - No:	TDS-CT50RS-L	Revision: A0
Remarks	single use only	
Specification		
Product name	CT50RS-L, 50ml Centrifuge tube	
Production conditions	Class 100,000 clean room production	
	No release agent is used	
Classification	CE/IVDD	
Dimensions		
Colour	Clear tube with Blue cap	
Graduation	Printing	
Nominal Capacity	50ml	
Material Characteristics	tube:polypropylene (PP)	
	cap:High Density Polyethylene(HDPE)	
Sterile	Yes	
Performance parameters	Storage Temperature:-80°C~121°C	
	Centrifugal force:10min 15,000xg	
Product Characteristics	biologically inert	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogen	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
	free of ATP	
Shelf - Life	60 months at controlled storage conditions	
Packaging	pcs per rack: 25	
	rack per case: 20	
Storage	dry condition:40%-60% rel. humidity	
Approved by PakGent quality department		