


Technical Datasheet

Article - No:	TDS-DI-L10	Revision: A0																										
Remarks	single use only																											
Specification																												
Part number	DI-L10																											
Product name	disposable inoculation loop, 10ul																											
Production conditions	No release agent is used																											
Classification	CE/IVDD																											
Dimensions																												
	<table border="1" style="width: 100%; border-collapse: collapse; text-align: center;"> <thead> <tr> <th style="font-size: small;">Description</th> <th style="font-size: small;">A</th> <th style="font-size: small;">B</th> <th style="font-size: small;">C</th> <th style="font-size: small;">D</th> <th style="font-size: small;">E</th> <th style="font-size: small;">F</th> <th style="font-size: small;">G</th> <th style="font-size: small;">H</th> </tr> </thead> <tbody> <tr> <td style="font-size: x-small;">Inoculation loop</td> <td style="font-size: x-small;">200</td> <td style="font-size: x-small;">130</td> <td style="font-size: x-small;">66.8</td> <td style="font-size: x-small;">0. D. 3. 6</td> <td style="font-size: x-small;">0. D. 6. 5</td> <td style="font-size: x-small;">0. D. 2. 7</td> <td style="font-size: x-small;">0. D. 3. 1</td> <td style="font-size: x-small;">0. D. 3. 8</td> </tr> <tr> <td style="font-size: x-small;">Inoculation loop</td> <td style="font-size: x-small;">198. 3</td> <td style="font-size: x-small;">130</td> <td style="font-size: x-small;">66. 8</td> <td style="font-size: x-small;">0. D. 0. 8</td> <td style="font-size: x-small;">0. D. 3</td> <td style="font-size: x-small;">0. D. 2. 7</td> <td style="font-size: x-small;">0. D. 3. 1</td> <td style="font-size: x-small;">0. D. 3. 8</td> </tr> </tbody> </table>		Description	A	B	C	D	E	F	G	H	Inoculation loop	200	130	66.8	0. D. 3. 6	0. D. 6. 5	0. D. 2. 7	0. D. 3. 1	0. D. 3. 8	Inoculation loop	198. 3	130	66. 8	0. D. 0. 8	0. D. 3	0. D. 2. 7	0. D. 3. 1
Description	A	B	C	D	E	F	G	H																				
Inoculation loop	200	130	66.8	0. D. 3. 6	0. D. 6. 5	0. D. 2. 7	0. D. 3. 1	0. D. 3. 8																				
Inoculation loop	198. 3	130	66. 8	0. D. 0. 8	0. D. 3	0. D. 2. 7	0. D. 3. 1	0. D. 3. 8																				
Colour	Orange																											
Material Characteristics	polypropylene (PP)																											
Shelf - Life	60 months at controlled storage conditions																											
Packaging	25pcs/bag, 40bags/case																											
Storage	dry condition: 40%-60% rel. humidity																											
Approved by PakGent quality department																												