


# Petri Dish


Petri dishes are one of the most important products in bacteriological laboratories. Manufactured from high-grade polystyrene, the petri dishes have exceptional optical clarity for microscopic analysis as well as heat resistance for use with hot agar. To enable parallel testing of a single sample using different agars, divided dishes with two or three or four compartments or quadrant grid are available.




Features	Description
<b>Material</b>	High-grade polystyrene
<b>Compatibility</b>	Compliant with automatic plate pourers
<b>Specification</b>	35 mm, 60 mm, 90 mm, 100 mm, 150 mm
<b>Appearance</b>	Divided dishes with two or three or four compartments, quadrant grid, flat, are available
<b>Sterilization method</b>	EO (Sterilization is optional)



Transparent and clean  
easy to observe under a microscope



Stable stacking



Double-layer protective packaging for  
easy use and transportation



Part NO.	Description	Pack
PD90	Petri Dish, 90*15mm, PS, whithout external grip	20 pcs/bag, 25 bags/case
PD90B	Petri Dish, 90*15mm, PS	10 pcs/bag, 50 bags/case
PD60	Petri Dish, 60*15mm, PS, whithout external grip	20 pcs/bag, 50 bags/case
PD35	Petri Dish, 35*15mm, PS, whithout external grip	10 pcs/bag, 230 bags/case
PD90-2	Petri Dish, 90*15mm, 2 compartments, PS, whithout external grip	20 pcs/bag, 25 bags/case
PD90-3	Petri Dish, 90*15mm, 3 compartments, PS, whithout external grip	20 pcs/bag, 25 bags/case
PD90-4	Petri Dish, 90*15mm, 4 compartments, PS, whithout external grip	20 pcs/bag, 25 bags/case
PD60-RGC	Rodac Dish, quadrant grid, 60mm, PS, concave, whithout external grip	20 pcs/bag, 50 bags/case
PD60-RGF	Rodac Dish, quadrant grid, 60mm, PS, flat, whithout external grip	20 pcs/bag, 50 bags/case
PD100	Petri Dish, 100*15mm, PS	10 pcs/bag, 32 bags/case
PD150	Petri Dish, 150*15mm, PS	10 pcs/bag, 16 bags/case

# PakGent®

Save every drop of precious reagent for you.

**PakGent Bioscience (Suzhou) Co., Ltd.**  
 Tel: +86 (0) 512 6670 0017  
 Address: Floor 1-3, Building 33, Alloy Industrial Park, Wuzhong District, Suzhou, China.

**PakGent Bioscience (Changshu) Co., Ltd.**  
 Tel: +86 (0) 512 5298 8330  
 Address: No.9 Lecheng Road, Bixi Street, Changshu, Jiangsu Province, China.

PakGent Bioscience (Shanghai) Co., Ltd.  
 Tel: +86 (0) 512 6670 0017  
 Address: Building C, No 888, Lake west road, Lin gang section, Free Trade Zone, Shanghai, China.