


Guidelines For Use


No. : PGLTT20231001E

Light Touch Pipette Tips

85% Less Force And Higher Accuracy




Without LTS
High operating force



With LTS
Low operating force

LTS pipette tips users can achieve reliable and consistent fitting effect without excessive force. LTS pipette tips has well-designed small sealing area and front stop point design. The front stop design allows you to know the exact position of the fitting.

Tip Ejection Forces




Please identify the LTS Pipette before purchasing

Multi-channel pipetting can be performed smoothly with the LTS tips system, ensuring absolutely consistent sample extraction across all channels. The front stop point design means that all channels are sealed perfectly at once, eliminating the need to shake or tighten the tips by hand. Significant time and money savings can be achieved during multi-channel pipetting, which ensures safe and reliable pipetting across all channels, consistent pipetting, and consistent, reproducible results.


20ul Fit LTS Pipette, sterile

Part NO.	With Filter	Low Retention	Sterile	Package Type	Pack
RT20R	-	-	✓	rack pack	96 pcs/rack, 50 racks/case
RT20R-L	-	✓	✓	rack pack	96 pcs/rack, 50 racks/case
RT20RF	✓	-	✓	rack pack	96 pcs/rack, 50 racks/case
RT20RF-L	✓	✓	✓	rack pack	96 pcs/rack, 50 racks/case
RT20-REF	-	-	✓	refill pack	96 pcs/rack, 10 racks/unit, 10 units/case
RT20-REF-L	-	✓	✓	refill pack	96 pcs/rack, 10 racks/unit, 10 units/case



200ul Fit LTS Pipette, sterile

Part NO.	With Filter	Low Retention	Sterile	Package Type	Pack
RT200R	-	-	✓	rack pack	96 pcs/rack, 50 racks/case
RT200R-L	-	✓	✓	rack pack	96 pcs/rack, 50 racks/case
RT200RF	✓	-	✓	rack pack	96 pcs/rack, 50 racks/case
RT200RF-L	✓	✓	✓	rack pack	96 pcs/rack, 50 racks/case
RT200B	-	-	-	bulk pack	1,000 pcs/bag, 20 bags/case
RT200B-L	-	✓	-	bulk pack	1,000 pcs/bag, 20 bags/case
RT200BF	✓	-	-	bulk pack	1,000 pcs/bag, 20 bags/case
RT200BF-L	✓	✓	-	bulk pack	1,000 pcs/bag, 20 bags/case
RT200-REF	-	-	✓	refill pack	96 pcs/rack, 10 racks/unit, 10 units/case
RT200-REF-L	-	✓	✓	refill pack	96 pcs/rack, 10 racks/unit, 10 units/case




300ul Fit LTS Pipette, sterile

Part NO.	With Filter	Low Retention	Sterile	Package Type	Pack
RT300R	-	-	✓	rack pack	96 pcs/rack, 50 racks/case
RT300R-L	-	✓	✓	rack pack	96 pcs/rack, 50 racks/case
RT300RF	✓	-	✓	rack pack	96 pcs/rack, 50 racks/case
RT300RF-L	✓	✓	✓	rack pack	96 pcs/rack, 50 racks/case
RT300B	-	-	-	bulk pack	1,000 pcs/bag, 20 bags/case
RT300B-L	-	✓	-	bulk pack	1,000 pcs/bag, 20 bags/case
RT300BF	✓	-	-	bulk pack	1,000 pcs/bag, 20 bags/case
RT300BF-L	✓	✓	-	bulk pack	1,000 pcs/bag, 20 bags/case
RT300-REF	-	-	✓	refill pack	96 pcs/rack, 8 racks/unit, 10 units/case
RT300-REF-L	-	✓	✓	refill pack	96 pcs/rack, 8 racks/unit, 10 units/case




1000ul Fit LTS Pipette, sterile

Part NO.	With Filter	Low Retention	Sterile	Package Type	Pack
RT1000R	-	-	✓	rack pack	96 pcs/rack, 50 racks/case
RT1000R-L	-	✓	✓	rack pack	96 pcs/rack, 50 racks/case
RT1000RF	✓	-	✓	rack pack	96 pcs/rack, 50 racks/case
RT1000RF-L	✓	✓	✓	rack pack	96 pcs/rack, 50 racks/case
RT1000B	-	-	-	bulk pack	1,000 pcs/bag, 10 bags/case
RT1000B-L	-	✓	-	bulk pack	1,000 pcs/bag, 10 bags/case
RT1000BF	✓	-	-	bulk pack	1,000 pcs/bag, 10 bags/case
RT1000BF-L	✓	✓	-	bulk pack	1,000 pcs/bag, 10 bags/case
RT1000-REF	-	-	✓	refill pack	96 pcs/rack, 8 racks/unit, 10 units/case
RT1000-REF-L	-	✓	✓	refill pack	96 pcs/rack, 8 racks/unit, 10 units/case



1200ul Fit LTS Pipette, sterile

Part NO.	With Filter	Low Retention	Sterile	Package Type	Pack
RT1200R	-	-	✓	rack pack	96 pcs/rack, 50 racks/case
RT1200R-L	-	✓	✓	rack pack	96 pcs/rack, 50 racks/case
RT1200RF	✓	-	✓	rack pack	96 pcs/rack, 50 racks/case
RT1200RF-L	✓	✓	✓	rack pack	96 pcs/rack, 50 racks/case
RT1200B	-	-	-	bulk pack	1,000 pcs/bag, 10 bags/case
RT1200B-L	-	✓	-	bulk pack	1,000 pcs/bag, 10 bags/case
RT1200BF	✓	-	-	bulk pack	1,000 pcs/bag, 10 bags/case
RT1200BF-L	✓	✓	-	bulk pack	1,000 pcs/bag, 10 bags/case



PakGent®

PakGent Bioscience (Suzhou) Co., Ltd.

Address: No. 9 Lecheng Road, Bixi Street, Changshu, Jiangsu Province, China

Tel: +86 (0) 512 6670 0017

Email: info@pakgentbio.com

Website: www.pakgentbio.com



Catalog