


Technical Datasheet

Article - No:	TDS-PD150	Revision: A0
Remarks	single use only	
Specification		
Product name	PD150, 150mm Petri dish, 15mm height	
Production conditions	Class 100,000 clean room production	
	No release agent is used	
Classification	CE/IVD	
Dimensions		
Colour	Clear	
Nominal Capacity	150mm diameter	
Material Characteristics	polystyrene (PS), premium surface	
Sterile	EO sterilization	
Product Characteristics	biologically inert	
	no heavy metals according to EC/94/62	
Purity Grade	free of RNase and DNase	
	free of pyrogen	
	free of PCR inhibitors	
	free of human DNA	
	free of endotoxins	
Shelf - Life	6 0 months at controlled storage conditions	
	free of ATP	
Packaging	pcs per bag: 10	
	bag per case: 16	
Storage	Storage Temperature: 20°C-25°C, humidity 40%-60%	
	Keep clean, ventilated and out of direct sunlight	
Approved by PakGent quality department		